

Portable Dental X-ray System



ABKAR MEDICAL
EQUIPMENTS L.L.C

Reliable Partner in Helathcare

+971 56 2351038
+971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira , Dubai , UAE



Technical Parameters

Tube Voltage: 65kV
Tube Current: 2mA
Time Exposure: 0.1-2.5s
High pressure generator: 40kHz DC
Rated Power: 400VA
Focal Spot To Skin Distance: 110mm
Focus Spot Size: 0.4mm(TOSHIBA)
Charger Input Voltage: AC100V-240V±10%
Radiation Leakage: <0.025mGy/h
Input Voltage: 25.2V DC
Supply Frequency: 50/60Hz+1%
Battery: 24.2V DC



DX 91AME