

# AED Automated External Defibrillator



Technology for Life



## HR-501 (AED)

Radian AED (Automated External Defibrillator) helps recover normal cardiac function upon acute heart disease, such as acute cardiac arrest, by providing cardiac function analysis and, if necessary, electroshock.

Heart Guardian™ as an AED, provides precise analysis of victims' body conditions and configures exact amount of shocks needed, promising safety and quicker responses.



### Warranty

5 years defibrillator  
4 years battery  
2 years pad



### Light Weight

1.9kg for easy portability



### Affordable



### Adult & Paediatric modes



**ABKAR MEDICAL**  
EQUIPMENTS L.L.C

*Reliable Partner in Helathcare*

+971 56 2351038  
+971 55 5546653

Abkar Medical Equipments L.L.C  
P.O. Box: 11037, Shop #4, Zarouni Building  
Opp. Azhar al Madina Hypermarket,  
Al Maktoum Hospital Road,  
Deira, Dubai, UAE



	HR-501
Product Type	AED Life Mode
ECG Analysis	8.5 Sec.
Hi-Cap Charging	8 Sec.
Shock Output	150J(Adults), 50J(Kids)
Data Storage	40 Min. (Embedded memory)
Data Transfer Protocol	Bluetooth 2.0
Control Buttons	Power, Shock
Battery Type	Li-Mn(12V DC 4.2A)
Shock Frequency	200 in Max
Pad Type	Electrode Pads, disposable and adhesive
Dimension	230(W) x 320(L) x 85(H)mm
Weight	1,9 kg
Warranty*	5 Yrs. (2 Yrs. for pads/ 4 Yrs. for Batteries)
Paid Repair Service	10 Yrs

The above HR-501 technical specification would be modified as its upgraded.

INM00025-SEPT13