

ZONCARE-M5

Full Digital Color Doppler Ultrasonic
Diagnostic System

First Choice of High Efficient Portable
Ultrasound System



 **ABKAR MEDICAL
EQUIPMENTS L.L.C**
Reliable Partner in Helathcare

 : +971 56 2351038
 : +971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira, Dubai, UAE

<http://www.zoncare.com>





Features

- 12 inch high resolution LED monitor with 90°rotation
- Two activated transducer connectors
- 2D,M, PW,CW,TDI, Elastic, AMM, CMM, 3D/4D
- THI, SRI, TSI,TCI , EFov ,Panoramic, HR flow, B-steer
- Standard Auto IMT measurement
- Support local and global amplification
- Support PW auto track and auto calculation
- Preposed USB ports, supporting multi-data transmission
- Color focus automatic tracking, Dual Live and Multi-synchronization technologies
- Multiple pregnancies measurement, Fetal growth curve, OB table are standard configured

Abundant clinical applications

ZONCARE-M5 has excellent performance and overall applications to meet user's different demands: Abdomen, OB/GYN, Cardiac, Urology, Vascular, Small parts, MSK, Nerve blocks etc.



Optional trolley and suitcase

ZONCARE-M5

Full Digital Color Doppler Ultrasonic Diagnostic System

Switching from a Black & White system to a highly efficient Color Doppler system has now been made easier and simpler than ever before, making it possible to reach beyond boundaries and transcend the ordinary when it comes to patient care. M5 is a brilliant choice for those who require high-quality image performance, ease of mobility, as well as affordability when it comes to advanced ultrasound imaging.

