



ABKAR MEDICAL
EQUIPMENTS L.L.C

Reliable Partner in Helathcare



Main feature:

- Monitor: 12.1" high solution LCD.
- Back-lit silica gel keyboard, photoelectric trackball and encoder, easy to operate.
- Abundant software : Abdomen, Obstetric, Gynecology, Cardiac, Urological and other part measurement.
- 28 grades depth adjustment under full visual field, maximum display depth is 240mm.
- Puncture guiding function, optional puncture probe, realizing ultrasound intervention requirement.
- 128 frames permanent storage, 256 frames cine loop.
- Various output ports can be support images and reports printing directly..
- DICOM3.0 network port, easily realizes network transmission function.

Full Digital Portable Ultrasound Scanner

Create the health with science and technology



Technical Specifications:

- Display depth (mm): 240 (max),16 levels adjustable.
- Monitor size:12 inch.
- Display mode: B、B+B、B+M、M、4B.
- Image gray scale:256 level.
- Cine loop: 809 frame (Max).
- Image storage: 32 frames.
- Scan angle: Adjustable.
- Scan depth: 40mm-240mm.
- Acoustic power: 2 steps.
- Dynamic range: 100dB-130dB.
- Probe elements : 96.
- Image flip: Up/down, left/right, black/white.
- Focus position: Adjustable.
- Focal space: 5 level.

Standard Configuration:

- Scanner Body: 1pc
- 3.5Mhz Convex Probe: 1pc
- Charger with a Set of wires: 1set
- Bottle of USG gel: 1pc
- User Manual: 1pc

Optional:

- Scanner Body: 1pc
- 3.5Mhz Convex Probe: 1pc
- Charger with a Set of wires: 1set
- Bottle of USG gel: 1pc
- User Manual: 1pc



ABKAR MEDICAL
EQUIPMENTS L.L.C

Reliable Partner in Helathcare

+971 56 2351038
+971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira , Dubai , UAE