



ZONCARE-P7

Full Digital Color Doppler Ultrasonic Diagnostic System

zoncare
MC CE ISO13485



 : +971 55 5546653
 : +971 56 2351038
 Abkar Medical Equipments L.L.C
 P.O. Box: 11037, Shop #4, Zarouni Building
 Opp. Azhar al Madina Hypermarket,
 Al Maktoum Hospital Road,
 Deira, Dubai, UAE

<http://www.zoncare.com>






ZONCARE-P7

Full Digital Color Doppler Ultrasonic Diagnostic System

P7 is the perfect match for high quality image performance, with added features such as Auto IMT, 4D, Elastography and Tissue Doppler Imaging and with the most competitive price in the industry.




Abundant Clinical Applications




Breast elastography




Vertebral artery PW



Cardiac M mode




Cardiac C mode




Liver 4B mode




Kidney C mode



Fetal Spine



fetal face--4D



Radial artery



Liver B mode




GB



CCA C mode

Features

- 15/19 inch high resolution LED monitor
- 8.4 inch touch screen
- Four activated transducer connectors
- Standard configured gel heating cup and keyboard
- 2D, M, PW, CW, TDI, Elastic, AMM, CMM, 3D/4D
- THI, SRI, TSI, TCI, EFov, Panoramic, HR flow, B-steer
- Support local and global amplification
- Support PW auto track and auto calculation
- Preposed USB ports, supporting multi-data transmission
- Multiple pregnancies measurement, Fetal growth curve, OB table are standard configured
- Color focus automatic tracking, Dual Live and Multi-synchronization technologies



Broadband Frequency Transducers

