



Binocular Biological
Ayaan
Microscope



microscopes

+971 56 2351038
+971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira, Dubai, UAE



**ABKAR MEDICAL
EQUIPMENTS L.L.C**

Reliable Partner in Healthcare

01 HEAD

Binocular straight tube (tilt at 45°)
360° rotation, To meet the omni-directional
observation requirements



02 NOSEPIECE

4 Holes
Equipped with objective 4X 10X 40X 100X

03 WORKING STAGE

The platform is equipped with
movable calipers and clamps



04 LED SPOTLIGHT

With single lens condenser,
Easy to adjust the light source

05 FOCUSING

Coaxial coarse & fine focusing adjustment



BIOLOGICAL MICROSCOPE

EXPLORE MICRO WORLD



BIOLOGICAL MICROSCOPE XSP-500SM

Head	Binocular straight tube (tilt at 45°)
Eyepiece	WF10X WF16X
Objective lens	Achromatic 4X 10X 40X 100X
Magnification	40X-1600X
Nosepiece	4 Holes
Stage	Double deck movable square platform
Light source	LED Light(Up and down the light source)
Collecting lens	Abbe condenser



