



Ayaan

microscope



+971 56 2351038
+971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira, Dubai, UAE



**ABKAR MEDICAL
EQUIPMENTS L.L.C**

Reliable Partner in Healthcare

XSP-500E

Specification	
Achromatic Objectives	4X、10X、40X (S)、100X (S、Oil)
Eyepiece	Wide Field Eyepiece : WF10X、WF16X
Viewing head	Compensation Free Binocular Head Inclined at 30°
Stage	Double Layer Mechanical Stage Size 125X115mm
Focusing	Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.002mm
Condenser	Abbe NA=1.25 with Iris Diaphragm & Filter
IIIumination	LED(1W) Lamp 220V or 110V ; (LED 3W lamp optional, Halogen bulb 6V/20W optional)

