



ABKAR MEDICAL
EQUIPMENTS L.L.C

Reliable Partner in Helathcare

F



Hygiene
Multifunctional
Excellent performance
Super valve



Dentist stool



Aluminium foot frame



Nylon silent wheel



PU upholstery



Green - 0#



Light Blue - 1#



Middle Blue - 3#



Dark Blue - 4#



Mango Yellow - 6#



Orange - 7#



Claret - 8#



Apple Green - 9#



Black - 10#



Red - 11#



Pink - 12#



Silver - 13#



Grey - 14#



Rose - 15#



Nacarat - 16#



Jewel Green - 18#



: +971 56 2351038

: +971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira, Dubai, UAE



High quality chair motor



Built out detachable suction filter



Multi-articulated Headrest





Assistant unit

- 3-way syringe(warm water)
- Regulable suction handle
- Assistant touch panel for chair movement and cuspidor
- Pre-poistion for built-in curing light

Product Description

Dental light



Bulb: 4*LED



Color temperature:4800-5500



Illumination:7000-38000Lux



Control:
Infrared sensor switch&gradually intensity adjustment





Dentist table

- Air-brak instrument arm
- Broad instrument table with plastic pad
- Five positions hygienic instrument holder
- Three way syringe with removable tip
- Main touch panel for chair movement and cuspidor
- LED X-ray viewer





Main switch for air water,electricity

Assures the safety after off duty,convenient for on/off and increase of service life



Water unit

- Rotabel & detachable ceramic spittoon
- Programmable cup filler and spittoon rinsing system
- Constant temperature warm system
- Pure water supply system
- Solenoid valve with built-in filter

Dentist stool

- Dentist stool (high quality)with backrest
- Gas spring moving system
- Moveable stool base



Standard Configurations:	
3-way syringe(cold/hot)	2 SETS
Powerful suction device	1 SET
Saliva ejector	1 SET
LED sensing light	1 SET
Warm water system	1 SET
Pure water supply system	1 SET
Rotatable ceramic spittoon	1 SET
DC 24V motor driven chair	1 SET
Multifunctional assisant device	1 SET
Dentist Stool	1 SET
LED x-ray viewer	1 SET



LED light

- Control: Infrared sensor switch&gradually intensity adjustment
- Color temperature: 4800-5500
- Build: 4 x LED

Dental led light, Infrared sensor switch&gradually



Assistant unit

- 3 way syringe(warm water)
- Regulable suction handle
- Assistant touch panel for chair movement and cuspidor
- Pre-position for built-in curing light

dental assistant unit, touch panel for chair movement and cuspidor. pre-polstion for built in curing light.



Dentist table

- Air-brak instrument arm
- Broad instrument table with plastic pad
- Five positions hygiene instrument holder
- Three way syringe with removable tip
- Main touch panel for chair movement and cuspidor
- LED X-ray viewer

Dentist table have three way syringe with removable tip, main touch panel for chair movement and cuspidor etc.



Saliva Ejector

- water and gas quickly through the strong and weak suction generator to form a negative pressure for the doctor's treatment when the saliva use.


Dental saliva ejector , water and gas quickly through the strong and weak suction generator to form a negative pressure for the doctor's treatment.



Water unit

- Rotabel & detachable ceramic spittoon
- Programmable cup filler and spittoon rinsing system
- Constant temperature warm system
- Pure water supply system
- Setenold valve with built-in filter

Water unit can rotabel and detachable ceramic spittoon



Emergency stop switch button

water tank gas water source switch.

Emergency stop switch button can be quickly paused.

