

KX5600

Full digital B mode ultrasound scanner

NEW




4.5Kg

Net weight of main unit is 4.5Kg, smart and light, convenient for carrying.



200min

Large-capacity lithium battery can drive the main unit to work for more than 200 minutes.



The hard disk can be extended to 4G to get 5000 frames permanent image storage.

KAI XIN

KUZHOU KAI XIN ELECTRONIC INSTRUMENT CO.,LTD.

Humanized design

- Rapidly one-key operation includes one-key storage, one-key browse and one-key print.
- Unique CPU temperature control function, when temperature is too high, it can automatically shut down.
- Rotatable LED Panel to make the users get the best image visual sense from different angles.
- One-key to light the keyboard, suitable for operation in dark room.
- Practical function knob design, use A key to realize quick adjustment in a variety of conditions.



Powerful data management system

- Optimization of the query image function, can realize 1-2 seconds response to your needs.
- Neoteric internal workstation multi-plan template, quickly finish data, report management function, realize process multi medical histories function.

Abundant software package

- Equipped with professional gravel package which can realize real-time positioning measurement.
- OB measurement and analysis: 15 kinds OB table and several GA measurement, AFI, Fetus growth curve and OB report.
- GYN measurement and analysis: general GYN software packages and GYN report.
- Urological measurement and analysis: residual urine, volume of prostate, Urological report.



High quality clinical image

- High resolution and full screen display, meet your vision needs.
- Adjustable LED brightness, convenient and practical.



Standard configuration:

- Main unit: 1 pc
- 3.5MHz multi-frequency convex probe: 1 pc

Optional configuration:

- 6.5MHz multi-frequency trans-vaginal probe
- 7.5MHz multi-frequency linear probe
- Lithium Battery
- Charger
- Trolley

KAI XIN
KUZHOU KAI XIN ELECTRONIC INSTRUMENT CO.,LTD.

AHMED SAVAD
Biomedical Engineer

+971 56 2351038

Shop #4, Zarouni Building,
Opp. Azhar Al Madeena Hyper Market,
Al Maktoum Hospital Road,
Deira, Dubai, U.A.E.
P.O.Box: 11037

056 2351038
055 5346659



ABKAR
MEDICAL
EQUIPMENTS L.L.C.

إبكر للمعدات الطبية ش.ذ.م.م

Reliable Partner in Healthcare

abkarmed@hotmail.com