


Suction Machine
7E-A



ABKAR MEDICAL
EQUIPMENTS L.L.C
Reliable Partner in Helathcare

☎ :+971 56 2351038
📞 :+971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira , Dubai , UAE



Suction Machine
7E-A



Vacuum pressure meter

Stepless pressure adjustment



Technical data.

Input power: 90VA

Power supply: DC 12V; 44VA

Max negative pressure value: $\geq 0.075\text{MPa}$

Pumping rate: $\geq 17\text{L/min}$ (7E-A)

Operation noise: $\leq 65\text{dB}$ (A)


Reservoir volume: 1000mL, 1piece



Concealed liquid storage bottle



Free-oil pump



Quiet treatment



Medical Grade



Stepless pressurize



Low noise



Concealed bottle



ABKAR MEDICAL
EQUIPMENTS L.L.C

Reliable Partner in Helathcare



:+ 971 56 2351038

:+971 55 5546653

Abkar Medical Equipments L.L.C
P.O. Box: 11037, Shop #4, Zarouni Building
Opp. Azhar al Madina Hypermarket,
Al Maktoum Hospital Road,
Deira, Dubai, UAE